



## Initialization

1. Check bottles



- Fill the water bottle with deionized water
- Check that the waste bottle is empty
- TIP: Adding a drop of soap to the waste bottle helps to keep it clean*

2. Switch instrument on

- Make sure the instrument is connected to your laptop via ethernet cable and ethernet-to-USB adapter
- The power switch is at the rear of the instrument

3. Copy template

- Create a copy of the measurement template and rename it. This will be your workspace where the measurements are saved.

4. Start AmphaSoft




- Double click the AmphaSoft 2.0 icon on the desktop

5. Open template


- In the AmphaSoft menu click on **Workspace > Load**. Navigate to the workspace, select the folder by clicking on it and click **Select Folder**.

6. Place chip




- Place the chip in the instrument using the 3 positioning pins (○) and gently push on the middle of the chip (▼) without touching the electrodes
- Carefully close the lid

7. Perform initial rinsing



- Click on **S\_1** in the left **Navigation** panel and make sure the correct chip number is configured in the measurement tab. If not, select the correct chip from the **Chip Id** dropdown menu.
- Click on **Instrument Control** in the **Navigation** panel
- Select the **Basic** tab and click on **Start Rinsing**  
The software will first connect to the instrument (green progress bar at the bottom) and then initiate the rinsing process
- After successful connection, a popup message appears. Make sure an empty container is placed below the sample holder and confirm by clicking **Continue Rinsing** twice
- TIP: If the water filter is dry or contains air bubbles, perform a few rinsing steps until the tubing is filled with deionized water*

8. Start measuring



- Return to the measurement view by selecting the first measurement (**S\_1**) in the **Navigation** panel on the left
- Enter the sample name and start with the sample preparation



## Sample Preparation

- 1 Buffer Preparation
- Take the buffer bottle out of the fridge. Check the buffer quality. It must be particle free and completely transparent.
  - Prepare a 50 ml buffer aliquot
  - Put the buffer bottle back into the fridge.
  - Let the aliquot equilibrate to room temperature
- TIP: The buffer recommendations list on [www.amphasys.com/download](http://www.amphasys.com/download) contains the recommended buffers by species*
- TIP: For hydrophobic pollen, add Tween20 (final concentration 0.05 %) to the buffer aliquot to facilitate suspension*


- 2 Pollen Collection
- PURE POLLEN**

  - Add the quantity of pollen for one measurement into an 1.5 ml Eppendorf tube
  - Recommended quantity:


**ANTHERS**

  - Add anthers into an 1.5 ml Eppendorf tube
  - Number of anthers:  
Wheat 3, Pepper 3, Potato 1, Brassicas 6, Cucumber 3, Melon 3, Watermelon 3

**WHOLE FLOWERS**


  - Add flowers into an 1.5 ml Eppendorf tube
  - Recommended quantities:  
Fennel: 1 umbellet  
Carrot: 1 umbellet
- 

Eppendorf Tube      Spatula
- NOTE: High cell concentration lead to a higher risk of chip clogging and inaccuracy of the data.*

- 3 Rehydration (if required)
- If the sample material is dehydrated, rehydrate for 30 min. Place the open Eppendorf tube with the sample in a rehydration box and close the box.
- TIP: The Quick Guide Measurement & Tips shows how to determine the optimum rehydration time*
- 

Rehydration Box

- 4 Buffer
- TIP: For anthers: If pollen cannot be released from the anther by shaking, slice them with a scalpel before adding the measurement buffer or squeeze them with tweezers after adding the buffer (step 5)*
- Add 1 ml measurement buffer

- 5 Extraction and Suspension
- Suspend the pollen by finger flicking or moving the tube over a tube rack (see image)
- TIP: For anthers: In case the pollen cannot be released from anthers by shaking, use tweezers to squeeze the anthers*
- 

move fast →

- 6 Filtration
- Use the recommended filter type (see buffer recommendations list on [www.amphasys.com/download](http://www.amphasys.com/download)). Filters must be clean and dry.

- 7 Dilution
- F chip: Add 1 ml measurement buffer
  - D chip: Add 1 ml measurement buffer
  - E chip: Add 2 ml measurement buffer

- 8 Equilibration
- Let the pollen equilibrate for 2 – 3 minutes in the buffer
- Exceptions: Sensitive pollen, such as wheat, brassica, artichoke. These samples must be measured immediately.
- TIP: The stability of the cells in the buffer can be tested using the experiments suggested in the Quick Guide Measurement & Tips, page 2.*

- 9 Mixing
- Slowly invert the FACS tube to distribute sedimented cells and immediately start measuring
- TIP: Seal the FACS tube with Parafilm before inverting*