

# APPLICATION NOTE

## YEAST CELL VIABILITY


Yeasts are single-celled fungus microorganisms that are widely used in biopharmaceutical production, food and beverage industry and as model organism in basic research. In order to control and improve biotechnological processes, monitoring of yeast viability and cell density are crucial. In order to achieve that, Amphasys developed a rapid and accurate method for yeast viability testing and viable and dead cell density determination using impedance flow cytometry (IFC).

### Straightforward 3-Step Workflow




### YEAST VIABILITY


Characterization of > 1000 cells / sec by impedance flow cytometry



### Yeast heat inactivation process



### Heat Inactivation




Amphasys AG  
Technopark Lucerne  
CH-6039 Root D4  
Switzerland



### Ampha X30 Impedance Flow Cytometer

- Rapid
- Accurate
- Reproducible
- Label-free
- Portable for on-site analysis



### Sample Preparation

- Dilution of yeast cells from solid agar plate or liquid culture in AF6 buffer (direct sampling from a bioreactor has also been successfully demonstrated)
- Filtration using a 10 um filter

### Instrument Setup

- Frequencies 2 / 12 MHz
- Triggering x positive
- Level ~ 0.01
- Modulation 5
- Amplification 6
- Demodulation 0
- Chip 30 um

### Analysis

- Signals at 12 MHz allow a discrimination between viable, apoptotic and dead yeast cells
- Identification of non-viable population using a heat-inactivated sample
- Application of gates to all samples
- Percentages of viable, apoptotic and dead cells are calculated automatically
- Possibility to create a measurement report (includes cell concentrations)